

Home Health Solutions

Product Catalogue

2017



Home Health Solutions



Table of contents

Blood Pressure Monitor	4 - 5
Nebulizer	6 - 7
Thermometer	8 - 9
Lifestyle	10 - 11
Baby Care	12

CHECKY^{ultra}

Blood Pressure Monitor - AUTOMATIC UPPER ARM

Order Code: M028701
Product No: 8698864007673

CHECKY ULTRA is checking your daily blood pressure accurately and thanks to the big bright LCD, all values can be easily read.

- ❶ Wide LCD screen
- ❷ Cuff Wrapping Detector
- ❸ Hypertension risk indicator
- ❹ Irregular Heartbeat Detector (IHB)
- ❺ Average of last 3 readings
- ❻ 2x120 memories



GUARDY

Blood Pressure Monitor - AUTOMATIC UPPER ARM

Order Code: M028721
Product No: 8698864007697

GUARDY can be used by 4 users and is monitoring the risk of Hypertension according to WHO.

- ❶ 120 memories with date and time for 4 users
- ❷ Irregular Heartbeat Detector (IHB)
- ❸ WHO risk category indicator
- ❹ Average of last 3 readings
- ❺ Auto power off



STILY

Blood Pressure Monitor - AUTOMATIC UPPER ARM

STILY helps the user to place the cuff correctly around the arm with the “Cuff Wrapping Detector” function. This function improves the fit of the cuff and therefore the measurements are more reliable and accurate.

Order Code: M029321
Product No: 8698864011809

- 1 Wide LCD screen
- 2 Cuff Wrapping Detector
- 3 Hypertension risk indicator
- 4 Irregular Heartbeat Detector (IHB)
- 5 Average of last 3 readings
- 6 2x120 memories



HANDY ultra

Blood Pressure Monitor - AUTOMATIC WRIST

HANDY ULTRA is a wrist watch device which can be easily hand-carried and is therefore ideally suitable for mobile use. Even in the dark, the attractive large backlit LCD display allows an easy reading of the accurately measured result.

Order Code: M028731
Product No: 8698864011748

- 1 Cuff Wrapping Detector
- 2 Hypertension risk indicator
- 3 Irregular Heartbeat Detector (IHB)
- 4 Average of last 3 readings
- 5 2x120 memories
- 6 Backlight LCD screen



BREATHY

Nebulizer - PISTON COMPRESSOR TYPE

Order Code: M029701
Product No: 8698864007628



The proprietary adjustable valve is able to deliver medications of different viscosity level according to every user's conditions without exchanging parts. The Valve Adjustable Technology (VAT) bottle allows users to adjust different levels of nebulization rate 0.15 (closed) / 0.5 (fully opened) ml/min at consistent particle size. Higher nebulization rate (fully opened) is for higher viscosity medications and higher breathing capacity user while lower nebulization rate with closed valve will be more appropriate for kids/infants with lower breathing capacity.

- 1 VA Technology: Patented valve adjustable technology
- 2 Particle size $< 3\mu\text{m}$
- 3 Effective cooling system for intensive use
- 4 Compact design
- 5 Low noise

VA
Technology

**Low Noise
Design**

**Consistent Fine
Particle Size
 $< 3\mu\text{m}$**



STEAMY

Nebulizer - COMPRESSOR TYPE

Order Code: M029721
Product No: 8698864007734



STEAMY is a simple and easy to use nebulizer with a powerful compressor. It's designed to deliver the prescribed medication solution to treat patient respiratory disorders, such as asthma, allergies and bronchitis. The nebulizer converts the medication solution into an aerosolized mist which is inhaled by the patient through the mouthpiece or masks.

- 1 Powerful compressor
- 2 Particle size $3\mu\text{m}$
- 3 Nebulizer kit
- 4 Compact design
- 5 Low noise



LUCKY

Nebulizer - COMPRESSOR TYPE

Order Code: M029723
Product No: 8698864012127



LUCKY has adorable, cute, child-friendly design. Operation of nebulizer is very easy with just one button. Suitable for adult and pediatric patients.

- 1 One button easy operation
- 2 Effective medication delivery
- 3 Adorable and friendly design
- 4 First choice, ideal for all children



SPEEDY

Digital Thermometer - NON-CONTACT INFRARED

Order Code: M029601
Product No: 8698864007611



SPEEDY non-contact temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.

- 1 1 second measurement (ASTM E1965-98:2009)
- 2 Body, object and room temperature
- 3 Fever alarm
- 4 9 memories
- 5 Backlight



SHINY

Digital Thermometer – PEN TYPE

SHINY is an anti-allergic thermometer. It does not include mercury and any metal components. You can measure the temperature in several parts on your body (Oral, axillar, rectal).

- 1 100% water resist
- 2 Fever alarm
- 3 Last measurement memory
- 4 Measurement start/finish alarm
- 5 Auto power OFF
- 6 4 colors



Order Code: M029621
Product No: 8698864010383

THERMY

Digital Thermometer – PEN TYPE

THERMY Digital Thermometer provides quick and highly accurate reading over the body temperature range. It is suitable for oral, rectal and axillary use.

- 1 $\pm 0.1^{\circ}\text{C}$ accuracy
- 2 Fast response
- 3 Memory function
- 4 Beeper function
- 5 Auto shut OFF



Order Code: M029602
Product No: 8698864011410

COSY

Order Code: M022170
Product No: 8698864010376

Anti-Decubitus AIR ALTERNATING PRESSURE MATTRESS SYSTEM

COSY air mattress is ideal for prevention and treatment of patients at low risk of pressure ulcers. It increases life quality of patients. Mattress is consist of a pressure adjustable pump for patients' weight.

- 1 Wavy pumping
- 2 Low noise and vibration
- 3 Adjustable pressure
- 4 Max. user capacity: 120 kg



SOOTHY

Reusable Pack - HOT / COLD

Order Code: M026003
Product No: 8698864011427

SOOTHY is used to treat pain by applying both hot and cold. The special covering provides comfortable treatment by preventing direct contact with the wounded area during both hot and cold applications.

- 1 It can be cooled or heated in microwave
- 2 Provides effective relief from muscular pain and helps reduce swelling
- 3 Great for treating a wide range of tissue injuries
- 4 Reusable



STICKY kit

Wound Care Kit

STICKY KIT provides two in one: practical skin cleansing and wound care with its alcohol cleaning stick and band-aid.

- 1 Cleaning stick containing 75% Ethyl alcohol
- 2 Band aid



Order Code: M026021
Product No: 8698864011793

STICKY

Cleaning Stick

STICKY is an effective alcohol cleaning stick ensures a quick & clean disinfection of the skin.

- 1 Cleaning stick containing 75% Ethyl alcohol



Order Code: M026020
Product No: 8698864011786

SMILKY

Breast Pump - ELECTRONIC

Order Code: M027501
Product No: 8698864007666

SMILKY electric breast pump is controlled by microcomputer chip. It has 2-phase suction which mimic sucking actions of infant. More milk for baby, more comfortable for Mum. The stimulating speed is fixed and stimulating suction divided into 9 adjustable levels. The expressing speed and suction are divided into 9 adjustable levels. It can mimic baby's sucking rhythm as well as sucking quantity. It can save or run the record of sucking speed and suction setting before hand.

- 1 Adjustable 9 speeds
- 2 Silicone cushion
- 3 BPA free
- 4 LCD screen
- 5 Memory function: Save the setting of last work under standby
- 6 Portable: Adapter & battery operated with lanyard and hook design
- 7 Quiet design

**BPA
Free**



Trimpeks' FIVE unique advantages at a glance!



01

Your 1-stop-shop

A complete product portfolio in both the HOME and the PROFESSIONAL segment allows you to experience **NEW HEALTH DIMENSIONS**. Several brands provide product and service solutions to meet the needs in the lifestyle, medical and professional fields.



02

Stock management and fast delivery

Our top modern warehouse and logistic center in Istanbul is handling all your orders professionally and quickly. Located in the center of EuroAsia, we guarantee fast delivery within 1-2 weeks after receipt of firm order. Where else do you get such leadtimes?



03

Small minimum order quantities

We offer very small MOQs without any risk and accordingly let dynamic and motivated customers focus on building-up a strong brand in the local market.



04

No Import Duty Tax

EU buyers benefit from orders accompanied by the «ATR.1 Free Circulation Document» which means no import duty tax for you.



05

Technical Service

Our very high-skilled Technical Service Dept. is responding to all issues promptly and with highest professionalism (24/7).

TRIMPEKS

new health dimensions

www.trimpeks.com

TRIMPEKS | Sultan Selim Mah. | Yunus Emre Cd.No:1/11 | 34415 Kağıthane | İSTANBUL
T (+90 212) 319 50 00 | F (+90 212) 319 50 50 | info@trimpeks.com