

PRODUCT CATALOGUE CATALOGUE PRODUITS



**Your “One-Stop Shop” for
Healthcare Requirements for**
✓Detection ✓Risk Assessment
✓ Diagnosis ✓Therapy



Trimpeks Healthcare has been established in 1990 in Istanbul/Turkey, by partners having the awareness of how valuable and important Health is, in a person's life.

Considering that a person should both manage to check himself and/or heal himself, at a daily routine or at the required times; we have built up a complete "Home Healthcare" range.

Considering as well that a person has at the required times, to be cared by the Health Professionals in Hospitals and Clinics; we have also built up a complete "basic" "Hospital Products" range.

In this way, we have become a **"1-Stop-Shop"** for healthcare requirements; which prestigiously serves over 50 countries around the world.

Our criteria for choosing the elements which build up our Product Portfolio is that:
"Each & every Trimpeks Product, should be of quality that can be used with our Family Members"

Since you have this Catalogue in your hand, we welcome you to be a Trimpeks Partner in your country.

Sincerely yours,

Your "TRIMPEKS TEAM"



INDEX

DIABETES CARE	4
TEMPERATURE MEASUREMENT	8
BLOOD PRESSURE MEASUREMENT	9
ACCESSORIES	11
STETHOSCOPES	12
PERCUSSION HAMMERS	13
NEBULIZERS	14
OXYGEN CONCENTRATORS	15
PULSE OXIMETERS	16
PAP DEVICES	17
CPAP MASKS	19
PRESSURE MATTRESS SYSTEMS	20
TENS & EMS ELECTRICAL STIMULATORS	21
PATIENT CARE	22

pM-100 *fasTtest*

Blood Glucose Monitor Device

Recognized for its ease of use and quick response time, the Plusmed fasTTest Blood Glucose Monitoring System is designed to process accurate results utilizing only a small amount of blood for in vitro diagnostics. The Plusmed fasTTest Blood Glucose Monitoring System is intended for use in the home and in professional settings to monitor whole blood glucose levels obtained from the fingertip, palm and forearm.



- 1 5 seconds measurement
- 2 0.5 µl blood sample
- 3 Auto- code system
- 4 4 reminder alarms
- 5 300 memories
- 6 Transfer to PC
- 7 Big LCD screen
- 8 Hypoglycemic & Hyperglycemic indicator alarms
- 9 7-14-30 days average results
- 10 "Hi" and "Lo" range settings
- 11 Alternative site testing

Order Code: M020751 - Product No: 5060186753614

pM-100 *fasTtest*

Blood Glucose Test Strips



The Plusmed fasTTest Blood Glucose Test Strips are used with the plusmed fasTTest Blood Glucose Meter for quantitatively measuring glucose (sugar) in whole blood obtained from the fingertip, palm and forearm. The Plusmed fasTTest Test Strips are for testing outside the body (in vitro diagnostic use).

Order Code: M020752 - Product No: 5060186753621

- 1 For testing glucose in capillary blood
- 2 Tiny blood sample size (0.5 µl)

pM-100 *fasTtest*

Blood Lancets



Plusmed fasTTest blood lancets are used with lancing device for taking blood sample. Plusmed fasTTest blood lancets are disposable and for single patient use.

Order Code: M021965 - Product No: 6957138609251

- 1 Painless application due to thinness
- 2 Needle gauge: 30G
- 3 More reliable due to gamma sterilization method

Recognized for its ease of use and quick response time, the Plusmed fasTest DailyCheck Blood Glucose Monitoring System is designed to process accurate results utilizing only a small amount of blood for in vitro diagnostics. The Plusmed fasTest DailyCheck Blood Glucose Monitoring System is intended for use in the home and in professional settings such as hospitals and clinics, to monitor whole blood glucose levels obtained from alternate sites in different positions such as sitting on a chair or at bedside. The display with contrast, large icons and large numbers enables people to read the test results easily.



Order Code: M020775 - Product No: 8698864016071

- 1 5 seconds measurement
- 2 1 µl blood sample
- 3 No coding required
- 4 Alarm for testing schedule
- 5 360 test memories with date&time
- 6 7,14 and 28 day averaging
- 7 mmol/L or mg/dL (switchable)

Order Code: M020776 - Product No: 8698864016088



pM-900 fasTest DiagnoPRO

Blood Glucose Monitoring System

Recognized for its ease of use and quick response time, Plusmed fasTest DiagnoPro Multi-Monitoring System is designed to process accurate results utilizing only a small amount of blood for in vitro diagnostics. Plusmed fasTest DiagnoPro Multi-Monitoring System is for quantitatively measuring Blood Glucose (sugar), Blood -Ketone and Blood Uric Acid in fresh capillary whole blood, neonatal blood, arterial and venous blood samples. It is only for use outside the body (For in vitro diagnostic use) and intended for self-testings of patients at home and professional use in hospitals and clinics to monitor blood glucose levels.



Order Code: M020776 - Product No: 8698864016088

- 1 Multi Function (Blood Glucose, Blood - Ketone* and Blood Uric Acid*)
- 2 GOD 0.7 µl Blood Sample
- 3 Average glucose in 7, 14, 30, 60 and 90-days
- 4 Results in 5 Seconds
- 5 Memory Up to 500 Tests
- 6 Strip Ejector
- 7 ISO 15197:2013 / EN ISO 15197:2015 Certified
- 8 PC Download via data transfer cable
- 9 Suitable for Neonatal use
- 10 6 months test strip expiry after vial's first opening
- 11 Measurement range indicator

Order Code: M020722 - Product No: 8698864014572



Fastest Hblyzer is intended to be used for quantitative determination of Hemoglobin (Hb) and calculated Hematocrit (Hct) in capillary and venous human whole blood. This system is intended for outside of the body (in vitro diagnostic use) by physical examination and by healthcare professional in a clinical situation such as anaemia.



- 1 Results in 15 seconds
- 2 10 µl small blood sample size
- 3 Memory stores 1000 reading with date and time
- 4 Switchable three measuring unit
- 5 Deluxe carry case
- 6 Large LCD
- 7 Auto power off
- 8 Uploading Data to PC via Plusmed Fasttest Hblyzer Software

Order Code: M020783 - Product No: 8698864014954



fastTest Hblyzer hemoglobin strip is used with fastTest Hblyzer hemoglobin meter as a Hemoglobin Analyzing System to quantitatively measure hemoglobin with capillary and venous whole blood. The Hemoglobin Analyzing System is intended for in vitro diagnostic use (for Medical Detection Mechanism) only (outside of the body).

Order Code: M020781 - Product No: 8698864014961

- 1 For testing hemoglobin with capillary and venous whole blood
- 2 Tiny blood sample volume (10 µL)
- 3 50 test strips

GLYCOHEMOGLOBIN ANALYZER

Unique Handheld HbA1c.
Analyzer for Better Diabetes Care

- 1 Portable Testing of HbA1c
- 2 Accurate Result with CV<3%
- 3 Easy and Fast Operation
- 4 Room Temperature Storage



Order Code: M020784 - Product No: 8698864015975



Order Code: M020785 - Product No: 6935784400033

FASTFINE

Insulin Pen Needles

Pen needles are used in conjunction with injection pens to deliver injectable medications into the body. Insulin pen needles come in a variety of lengths and widths to suit all body types. From 8 mm to 4 mm, the needle length you choose is likely to be dictated by your size, children being likely to benefit from the shorter 4 mm size. Each pen needle is embedded in a plastic hub that can attach to the pen.



- 1 Sterile
- 2 Fits all insulin pens
- 3 Non toxic, non pyrogenic
- 4 Silicon coated
- 5 Available in 4 different sizes
 - 32 G (0.23 X 4 mm)
 - 32 G (0.23 X 5 mm)
 - 31 G (0.25 X 6 mm)
 - 31 G (0.25 X 8 mm)

4MM-32G-100 - Order Code: M026620 - Product No: 6972857290863

5MM-32G-100 - Order Code: M026621 - Product No: 6972857291723

6MM-32G-100 - Order Code: M026622 - Product No: 6972857290856

8MM-31G-100 - Order Code: M026623 - Product No: 6972857291716

MT-811

Pen Type - Digital Thermometer

This thermometer can measure the temperature in various locations on the body which maintain a fairly stable temperature (oral, axillary, rectally). Display readings to great precision (0.1 °C or 0.2 °F).

- 1 ± 0.1°C Accuracy
- 2 Fast response
- 3 Memory function
- 4 Auto shut-off

Order Code: M029430 - **Product No:** 6973709231058



UltraScan II

Infrared Thermometer Ear & Forehead Scan

- 1 Measured on forehead or ear
- 2 Fast reading
- 3 10 reading memory
- 4 3-color Backlight
- 5 Result in 1 second
- 6 Automatic power-off

Order Code: M020438 - **Product No:** 6970392214481



Recognized for its ease of use and quick response time, the Plusmed fastTest DailyCheck Blood Glucose Monitoring System is designed to process accurate results utilizing only a small amount of blood for in vitro diagnostics. The Plusmed fastTest DailyCheck Blood Glucose Monitoring System is intended for use in the home and in professional settings such as hospitals and clinics, to monitor whole blood glucose levels obtained from alternate sites in different positions such as sitting on a chair or at bedside. The display with contrast, large icons and large numbers enables people to read the test results easily.

- 1 Accurate measurements
- 2 Big Screen with easy to read Big Numbers and Backlight
- 3 WHO Hypertension Risk Indicator
- 4 Irregular Heartbeat detector (IHB)
- 5 Average of the last 3 measurements
- 6 2 x 120 memories with date & time indication for two users
- 7 Cuff Tightness indicator
- 8 Body Movement indicator
- 9 Low Battery indicator
- 10 Auto power off
- 11 Type-C USB port
- 12 22-42 cm Cuff
- 13 4 x "AAA" Alkaline Battery
- 14 Carry Pouch



Order Code: M028757 - **Product No:** 6948158401154

CONNECT

Arm-Type Fully Automatic Digital Blood Pressure Monitor

- 1 Accurate Measurement
- 2 Quiet Inflation Method
- 3 Automatic Power Off
- 4 Irregular Heartbeat Indicator
- 5 Last 3 Tests Average
- 6 Static Pressure Measurement
- 7 Bluetooth connection
- 8 2x150 Memories in Two Groups with Date and Time
- 9 Cuff Circumference Approx.135 (W) x 485(L) mm



Order Code: M028743 - **Product No:** 6970392213408

pM-8197W

Wrist-Type Fully Automatic Digital Blood Pressure Monitor

- 1 Quiet Inflation Method
- 2 Average of Last 3 Measurements
- 3 Accurate Measurement
- 4 Irregular Heartbeat Detection
- 5 Quick Result
- 6 Automatic Power-off
- 7 Wrist Circumference: 13,5-21,5 cm (5,3-8,5")



Product No: 8698864016484

PERFECT

Professional Aneroid Sphygmomanometer

- 1 Chrome plated metal case
- 2 48 mm diameter of dial
- 3 Stainless steel with chrome plated
- 4 0-300 mmHg scale
- 5 ± 3 mmHg measurement accuracy
- 6 Green Cotton Adult Hook Cuff
- 7 Double tube system
- 8 Storage case

Order Code: M021629 - **Product No:** 6970384252132



Palm Sphyg I

Aneroid Sphygmomanometer

- 1 ABS plastic upper case
- 2 60.2 mm diameter of dial
- 3 Right or left hand usage
- 4 Large zinc-alloy handle with black powder-coating
- 5 Deflation valve: "Push to button" for rapid deflation, "screwing" for precise deflation
- 6 Nylon D-Ring Adult Cuff
- 7 Silicone shockproof ring
- 8 Storage case

Order Code: M021630 - **Product No:** 6970384252187



pM-A01S

Aneroid Sphygmomanometer

Blood pressure measurements with the PM-A01S are done using the cuff / stethoscope listening method. It provides more reliable measuring.

- 1 300 mmHg non-stop pin manometer
- 2 D-ring nylon cuff
- 3 Heavy duty metal gauge
- 4 Lightweight exterior single head stethoscope
- 5 Chrome plated brass air relief valve

Order Code: M029350 - **Product No:** 6975531370097



METAL HOOK CUFF

- 1 Chrome plated metal case
- 2 48 mm diameter of dial
- 3 Stainless steel handle with chrome plated
- 4 0-300 mmHg scale
- 5 ± 3 mmHg measurement accuracy
- 6 Green Cotton Adult Hook Cuff
- 7 Double tube system
- 8 Storage case



MANOMETER

- 1 Compatible with Aneroid
- 2 Sphygmomanometers



BULB

- 1 Compatible with Aneroid
- 2 Sphygmomanometers



VELCRO CUFFS

Specification	Size (mm)	Arm Range (mm)
New born	185x55	----
Infant	255x75	130~190
Child	345x110	184~267
Adult	540x145	254~406
Large adult (Obese)	620x175	333~510
Thigh	780x220	406~660



D-RING CUFFS

Specification	Size (mm)
Adult	500x140
Large adult (Obese)	600x170



pM-ST71

Stethoscope - Single Head

- 1 Comfortable fit
- 2 High Quality
- 3 Better healthcare
- 4 Good acoustic response
- 5 Ear tips with soft cone-shaped bell
- 6 Aluminum alloy
- 7 Aluminum alloy binaural with chromium plate
- 8 Chest piece diameter: 46 mm
- 9 Diaphragm diameter: 42.8 mm

Order Code: M020180 - **Product No:** 8698864015524



pM-ST72

Stethoscope - Dual Head

- 1 Comfortable fit
- 2 High Quality
- 3 Better healthcare
- 4 Good acoustic response
- 5 Bell with non-chilling ring
- 6 Ear tips with soft cone-shaped bell
- 7 Opposite direction, two feature function
- 8 Aluminum alloy
- 9 Aluminum alloy binaural with chromium plate
- 10 Chest piece diameter: 46 mm
- 11 Diaphragm diameter: 42.8 mm

Order Code: M020195 - **Product No:** 8698864015531



CardioSTETH I

Stainless Steel Cardiology Stethoscope

- 1 Silver dual head chestpiece for cardiology diagnosis
- 2 Stainless steel chestpiece, excellent sound conduction
- 3 The membrane is sensitive for sound collecting
- 4 A" shaped thick tube with excellent resistance to ambient
- 5 interference ensuring best sound conduction
- 6 Comes complete with accessories

Order Code: M020190 - **Product No:** 6970384252194



TaylHM

Taylor Type Percussion Hammer

- ① Weight: 56 g
- ② Length: 188 mm
- ③ Width: 21 mm
- ④ Zinc alloy handle, PVC head



Order Code: M030035 - Product No: 6970384252200

BuckHM

Buck Type Percussion Hammer

- ① Weight: 90 g - Length: 186 mm - Width: 56 mm
- ② Zinc alloy handle, PVC head
- ③ Diagnostic pin can be screwed off conveniently in the head, while diagnostic brush can be taken out in the end



Order Code: M030036 - Product No: 6970384252217

BabiHM

Babinski Type Percussion Hammer

- ① Weight: 97 g
- ② Length: 230 mm
- ③ Width: 44 mm
- ④ Zinc alloy handle, PVC head



Order Code: M030037 - Product No: 6970384252224

RespiAid SMART

Piston Compressor Nebulizer

- 1 Piston Compressor
- 2 Ample Ventilation
- 3 Safety Fuse
- 4 2-Piece Medicine Bottle
- 5 Low Noise Level - Less than 60 dBA
- 6 Wall-Mountable Design
- 7 Shock-Resistant Case with Handle
- 8 3-Ball Bearing for Smoother Operation

Order Code: M029562 - Product No: 8906099802520



RespiAid Mesh-01

Mesh Nebulizer

Plusmed RespiAid Mesh-01 is intended for the treatment of respiratory diseases such as asthma, allergies, etc. Thanks to the self-excited oscillation circuit, the ceramic sheets contained in nebulizing head together with the metal mesh are driven successively to oscillate at a high frequency, which makes the liquid flowing through the micromesh of the metal mesh and being ejected and nebulized to form an inhalant aerosol.

- 1 Replaceable medicine cup
- 2 Quiet design
- 3 Dual power supply (2 x 1.5 V AA batteries or Micro USB 5 V)
- 4 Particle size $\leq 5\mu$
- 5 Nebulization rate ≥ 0.20 ml/min
- 6 Capacity of medicine cup ≤ 8 ml
- 7 Accessories: Medicine Cup | Mouthpiece | Adult Mask | Child Mask | USB Cable |
- 8 Storage Bag | Battery

Order Code: M0295572 - Product No: 8698864016378



RespiAid P2

Portable Oxygen Concentrator

Plusmed RespiAid P2 Portable Oxygen Concentrator is used on a prescriptive basis by patients requiring supplemental oxygen. It supplies a high concentration of oxygen to the patients.

It supplies a high concentration of oxygen and is used with a nasal cannula to channel oxygen from the concentrator to the patient. The Plusmed RespiAid P2 is small, portable and may be used in home, institution and various mobile environments. Plusmed RespiAid P2 Portable Oxygen Concentrator conforms to all applicable FAA acceptance criteria for Portable Oxygen Concentrator carriage and use on board aircraft.

- ① Lightweight 1.98kg
- ② Mini Size: 221*85*160mm, Very small size, easy to carry
- ③ 5 Setting Flows: From 1-5 setting flows, can meet different oxygen demand
- ④ Pulse Oxygen Supply: Plusmed RespiAid P2 is equivalent to 5 liters
- ⑤ High Efficiency Molecular Sieve: Efficient adsorption capacity,
- ⑥ can purify more pure
- ⑦ Oxygen Concentration: 90% -3% / +6%
- ⑧ Powerful Compressor: Powerful compressor can provide stable output
- ⑨ Triple Filtration: Effectively remove air pollen, large particles, dust and other impurities
- ⑩ Battery: 8 cell standard Battery Up to 4 hours



Order Code: M023033 - Product No: 8698864015685

pM-KN04

Oxygen Concentrator

PM-KN04 Oxygen Concentrator consisted of air compressor, adsorption tower of molecular sieve, sensor and flow meter. With the molecular sieve to be absorbent, the machine adopts PSA (pressure swing adsorption) method to directly concentrate medical oxygen with high concentration from the air.

It is small and portable. The machine takes the advantages of small in volume, light in weight, convenient in moving (as it has turning truckle), stable in performance, high in safety, easy in operating, low in noise, safe in working, which complies with requirements of medical device. The machine can also realize the measurement and display of SpO₂, and pulse rate by connecting accessories. The machine has a wide application, including family, clinic, hospital, oxygen bar, plateau outpost, gym center, etc.

- ① With working indicator
- ② Built-in filters to pure the gas
- ③ Oxygen concentration status indicator
- ④ Function of accumulating time
- ⑤ Function of timing shutdown
- ⑥ Alarm for power failure
- ⑦ Alarm for high/low pressure protection
- ⑧ Alarm for low flowrate
- ⑨ Alarm for high temperature
- ⑩ Measure SpO₂ and pulse rate
- ⑪ Internal nebulizer



Order Code: M023027 - Product No: 8698864015692

OxyCheck

Pulse Oximeter

OxyCheck Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO₂) and the pulse rate of adult and pediatric patients in home and hospital environments (including clinical use in internist/-surgery, anesthesia, intensive care etc.).

- 1 Big display
- 2 Measures blood oxygen saturation level and pulse rate
- 3 SpO₂ & PR waveforms
- 4 Perfusion Index (PI)
- 5 Easy to Operate
- 6 High Accuracy
- 7 Low battery indicator
- 8 Automatic power off within 5 seconds
- 9 Screen displayed in 4 directions



Order Code: M061607 - **Product No:** 8698864015562

pM-AF03

Hydrophobic Suction Filter

- 1 Size: H:54mm; w: 64mm
- 2 Housing:ABS
- 3 Filter: Hydrophobic PTFE Media Disc
- 4 BFE: 99.999%
- 5 Pore Size:0.2µm
- 6 Weight:18.2g
- 7 Connector:8mm
- 8 Stepped Hose Barbs
- 9 Max Operating Pressure : 20psi
- 10 Effective Filtration Area: 24.6cm²
- 11 Max Operating Temp:60°C



Order Code: M024911 - **Product No:** 8698864016163

pM-AF05

Hydrophobic Suction Filter

- 1 Size: H:54mm;w:58mm
- 2 Housing:PPFilter:
- 3 Air filtration membrane
- 4 BFE:99.999%
- 5 Pore Size:0.3µm
- 6 Weight:19.8g
- 7 Connector:6mm-8mmStepped Hose
- 8 Barbs
- 9 Max Operating Pressure:20psi
- 10 Effective Filtration Area: 15.5cm²
- 11 Max Operating Temp: 60°C



Order Code: M024912 - **Product No:** 8698864016170

RespiAid

CPAP System

RespiAid CPAP is equipped with required alarm and warning systems to offer simplicity, safety and comfort to the user to provide continuous treatment.

- | | |
|------------------------------|---|
| 1 Modes: CPAP | 12 Humidifier: |
| 2 Pressure Range: 4 – 20 hpa | Easy clean I Easy to fill water chamber I |
| 3 Ramp time: 0-60 minutes | Metal Material I Fast Heating I Good humidification |
| 4 LCD Display | 13 Noise Level: <30 dB |
| 5 Powerful Analysis Software | 14 Language: Multiple |
| 6 SD card | 15 Weight: |
| 7 Auto Altitude Compensation | Without humidifier : <1 kg |
| 8 Auto Leak Compensation | Humidifier included : <1.6 kg |
| 9 Auto on/off | 16 Dimensions: |
| 10 High performance motor | 140 mm × 156mm × 94 mm (without Humidifier) |
| 11 Humidification: 0-5 | 262 mm × 156mm × 94mm (with Humidifier) |



Order Code: M022507 - Product No: 8698864015647

RespiAid

AutoCPAP System

Designed to provide the user with a comfortable sleep, RespiAid Auto CPAP offers all kinds of details for tracking the treatment parameters, and as an extremely technological product, it is equipped with alarm and warning systems and provides simplicity, safety and comfort to the user.

- | | |
|------------------------------|---|
| 1 Modes: CPAP, APAP | 12 Humidifier: |
| 2 Pressure Range: 4 – 20 hpa | Easy clean I Easy to fill water chamber I |
| 3 Ramp time: 0-60 minutes | Metal Material I Fast Heating I Good humidification |
| 4 LCD Display | 13 Noise Level: <30 dB |
| 5 Powerful Analysis Software | 14 Language: Multiple |
| 6 SD card | 15 Weight: |
| 7 Auto Altitude Compensation | Without humidifier : <1 kg |
| 8 Auto Leak Compensation | Humidifier included : <1.6 kg |
| 9 Auto on/off | 16 Dimensions: |
| 10 High performance motor | 140 mm × 156mm × 94 mm (without Humidifier) |
| 11 Humidification: 0-5 | 262 mm × 156mm × 94mm (with Humidifier) |



Order Code: M022508 - Product No: 8698864015654

Plusmed RespiAid BPAP-30T “S BPAP” can be comfortably used during central apnea besides supporting spontaneous breathing. While providing respiration with minimum effort thanks to the improved breathing compatibility technology without tiring the user. Furthermore, it supports the healthcare professionals with detailed analysis and an integrated reporting program.

- | | |
|-------------------------------------|---|
| 1 Modes: CPAP, APAP, BPAP-S, BPAP T | 13 Humidifier: |
| 2 Pressure Range: 4 – 30 hpa | Easy clean I Easy to fill water chamber I |
| 3 Ramp time: 0-60 minutes | Metal Material I Fast Heating I Good humidification |
| 4 LCD Display | 14 Noise Level: <30 dB |
| 5 Powerful Analysis Software | 15 Language: Multiple |
| 6 SD card | 16 Weight: |
| 7 EP Flex | Without humidifier : <1 kg |
| 8 Auto Altitude Compensation | Humidifier included : <1.6 kg |
| 9 Auto Leak Compensation | 17 Dimensions: |
| 10 Auto on/off | 140 mm × 156mm × 94 mm (without Humidifier) |
| 11 High performance motor | 262 mm × 156mm × 94mm (with Humidifier) |
| 12 Humidification: 0-5 | |

Order Code: M022509 - Product No: 8698864015661



Plusmed RespiAid BPAP-30T Pro “ST BPAP” can be comfortably used during central apnea besides supporting spontaneous breathing. While providing respiration with minimum effort thanks to the improved breathing compatibility technology without tiring the user. Furthermore, it supports the healthcare professionals with detailed analysis and an integrated reporting program.

- | | |
|--|---|
| 1 Modes: CPAP, APAP, BPAP-S, BPAP T, BPAP-ST | 13 Humidifier: |
| 2 Pressure Range: 4 – 30 hpa | Easy clean I Easy to fill water chamber I |
| 3 Ramp time: 0-60 minutes | Metal Material I Fast Heating I Good humidification |
| 4 LCD Display | 14 Noise Level: <30 dB |
| 5 Powerful Analysis Software | 15 Language: Multiple |
| 6 SD card | 16 Weight: |
| 7 EP Flex | Without humidifier : <1 kg |
| 8 Auto Altitude Compensation | Humidifier included : <1.6 kg |
| 9 Auto Leak Compensation | 17 Dimensions: |
| 10 Auto on/off | 140 mm × 156mm × 94 mm (without Humidifier) |
| 11 High performance motor | 262 mm × 156mm × 94mm (with Humidifier) |
| 12 Humidification: 0-5 | |

Order Code: M022510 - Product No: 8698864015678



pM-FM

Full Face CPAP Mask

pM-FM Face Mask is comfortable flexible and quiet as it is designed to minimize the contact with the skin, which provides a comfortable feeling to the user during the treatment.

- 1 Anti-asphyxia valve
- 2 Adjustable silicone forehead support
- 3 Two sided ventilation
- 4 Oxygen entry port
- 5 Latex free
- 6 Adjustable Clips
- 7 Revolver head elbow

MEDIUM | Order Code: M027181 - Product No: 8698864011687

LARGE | Order Code: M027183 - Product No: 8698864011687



pM-NM

Nasal CPAP Mask

pM-NM Nasal Mask is designed for a comfortable, quiet and flexible use. Due to the minimized contact with the skin, the user has a comfortable feeling during the treatment.

- 1 Adjustable silicone forehead support
- 2 Ventilation
- 3 Oxygen entry port
- 4 Does not Contain Latex
- 5 Unbreakable clips
- 6 Revolver head elbow

MEDIUM | Order Code: M027180 - Product No: 8698864011670

LARGE | Order Code: M027182 - Product No: 8698864011670



pM-AD01

Tubular Pressure Mattress System

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. For emergent CPR, please pull the CPR yellow strip outside quickly to reduce the air fast.

- 1 Air ventilation function
- 2 Tubes on A-B lines provide alternating inflating
- 3 CPR - Fast Air Release Valve
- 4 Made of heavy duty, premium PVC
- 5 Easy to use
- 6 Adjustable
- 7 Low noise, low vibration, low temperature
- 8 Easy fixing on bed frame with hooks

Order Code: M022143 - **Product No:** 6975397940014



pM-2012 PLUS

Bubble Pressure Mattress System

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. It also provides individual home care setting and long term care of whom suffering from pressure ulcer and pain management as prescribed by a physician.

- 1 Easy to use
- 2 Adjustable
- 3 Low noise
- 4 Low vibration, low temperature

Order Code: M022185 - **Product No:** 6975397940007



pM-TE01

TENS & EMS Electrical Stimulator

Plusmed TENS+EMS Dual Therapy System is a dual channel digital TENS/NMES electrical stimulator for active treatment application. Transcutaneous Electrical Nerve Stimulator (TENS) can be used for symptomatic relief of chronic intractable pain. Neuromuscular electrical stimulation (NMES) can be used for relax muscle spasms, prevent or retard disuse atrophy, increase local blood circulation, re- educate muscles or maintain / increase the range of motion.

- 1 Digital technology
- 2 2 channels
- 3 8 different programs
- 4 Low battery indicator
- 5 12-step-timer
- 6 Adjustable intensity and treatment time

Order Code: M026502 - **Product No:** 8698864013544



AS9020

TENS & EMS Electrical Stimulator

- 1 TENS - EMS - MASSAGE functions
- 2 2 output with 24 modes
- 3 20 level intensity
- 4 20 minutes default timer
- 5 Adjustable from 10-60 minutes
- 6 Screen with backlight
- 7 Rechargeable lithium battery
- 8 Battery indicator (shows when the battery needs to be recharged)
- 9 Type-C charging interface

Order Code: m026503 - **Product No:** 6970890320783



pM-CH320

Reusable Hot / Cold Pack

The use of heat and cold is a standard treatment for a variety of injuries. The type and age of the injury or condition indicates whether hot or cold is appropriate. The use of hot and cold packs is one of the simplest methods of applying heat or cold to the affected area. Cold pack decrease the bruising associated with sprains, as well as the swelling and help to relieve the pain associated with these conditions by numbing the area. Hot pack relaxes muscles and improves blood circulation and produces a soothing effect on the body, further aiding in the healing process.

- 1 Comfortable and convenient application
- 2 Provides effective relief from muscular pain
- 3 Helps reduce swelling
- 4 Great for treating a wide range of tissue injuries
- 5 Reusable pack is perfect for continuous treatment at home



Order Code: M022141 - Product No: 8698864007123



Premium Condoms


- 1 Latex Rubber Condoms
- 2 Unit pack size: 1 Condom
- 3 Nominal width in mm: 52±2mm
- 4 Quantity: 100 Condoms
- 5 No method of contraception can provide 100% protection against pregnancy or the transmission of HIV or other STD's.
- 6 Made of natural rubber latex.
- 7 Do not use if package is damage.

Order Code: M030718 - Product No: 6970533250491



COLOSTOMY - ILEOSTOMY - UROSTOMY SYSTEM

- 1 One-Piece Colostomy Pouch (Closed or Drainable)
- 2 One-Piece Urostomy Pouch (With Tap)
- 3 Two-Piece Colostomy Pouch (Closed or Drainable)
- 4 Baseplate for Two-Piece Colostomy System
- 5 Baseplate for Two-Piece Urostomy System





INDEX

[illegible]

TRULY

Digital Thermometer - INFRARED NON-CONTACT

It measures the ambient and body temperature without contact with colored warning symbols.

- ① 0,5 second measurement
- ② LCD display with backlight
- ③ Fever alarm
- ④ 32 memories
- ⑤ °C or °F switchable
- ⑥ Automatic shut-off
- ⑦ Color Temperature Warner



Order Code: M029631 - **Product No:** 8698864013742

GUARDY PRO

Digital Blood Pressure Monitor

Guardy Pro can be used by 2 users and is monitoring the risk of Hypertension according to WHO.

- ① Accurate measurements
- ② Big Screen with easy to read Big Numbers and Backlight
- ③ WHO Hypertension Risk Indicator
- ④ Irregular Heartbeat detector (IHB)
- ⑤ Average of the last 3 measurements
- ⑥ 2 x 120 memories with date & time indication for two users
- ⑦ Cuff Tightness indicator
- ⑧ Body Movement indicator
- ⑨ Low Battery indicator
- ⑩ Auto power off
- ⑪ Type-C USB port
- ⑫ 22-42 cm Cuff
- ⑬ 4 x "AAA" Alkaline Battery
- ⑭ Carry Pouch



Order Code: M028758 - **Product No:** 6948158401147

STEAMY

Nebulizer - COMPRESSOR TYPE

STEAMY is a user friendly nebulizer with a powerful compressor. It is designed to deliver the prescribed medication solution to treat patient respiratory disorders, such as asthma, allergies and bronchitis. The nebulizer converts the medication solution into an aerosolized mist which is inhaled by the patient through the mouthpiece or masks.

- 1 Average Nebulization Rate > 0.32ml/min
- 2 Mode of Operation: 24 hours continuous
- 3 Flow Rate: 5-8 LPM
- 4 MMAD: 3 µm

Order Code: M029724 - Product No: 869886400734



FONDLY

Electronic Baby Scale

- 1 4 high precision strain gauge sensor
- 2 Auto on and auto off
- 3 Tare function
- 4 Auto lock/hold function
- 5 Low battery and overload indication
- 6 Unit conversion switch
- 7 With Height Measurement Tape

Order Code: M024324 - Product No: 869886401573



BRP-306

Electric Breast Pump

- 1 Wearable type, Invisible, Handsfree Use
- 2 3 modes Massage Stimulation - Mixed Suction - Expression
- 3 9 selectable levels for each mode
- 4 Breast shield Size 24mm.
- 5 Built-in lithium battery with Type C cable included.
- 6 Display LCD and touch buttons to choose functions
- 7 Quiet design | Eco-sustainable packaging
- 8 1 x 17mm Flange Insert, 1 x 21mm Flange Insert

Order Code: M027503 - Product No: 6976313750236



WF401

Oral Irrigator

Power Supply	5V --- 1A
Water Tank	300ml
Power Consumption	5W
Battery Capacity	1800mAh
Charging Time	More than 4 hours

Order Code: M025108 - **Product No:** 6942284602010



Y4

Massage Gun

Size	154x121x57.3mm
Net weight	About 457g
Battery capacity	2x1800mAh
Charging interface	Type-C
Charging time	About 3.5 hours
Endurance time	About 6 hours
Maximum (blocking) current	7A
Gear/speed	6 Gear/2000 ~ 3000 RPM
Stroke	6mm
Maximum Power	About 55 W



Order Code: M025107 - **Product No:** 8698864016927



MEDIZINTECHNIK

INDEX

[illegible]

GENTLE

Wrist-Type Fully Automatic Digital Blood Pressure Monitor

- 1 Average of Last 3 Measurements
- 2 Accurate Measurement
- 3 Irregular Heartbeat Detection
- 4 Quick Result
- 5 Quiet Inflation Method
- 6 Automatic Power-off
- 7 Wrist Circumference: 13,5-21,5 cm (5,3-8,5")

Order Code: M028739 - Product No: 6942284616970



INTELLIGENT^{PRO}

Talking Digital Blood Pressure Monitor

- | | |
|--|---|
| <ol style="list-style-type: none"> 1 Accurate measurements 2 Big Screen with easy to read Big Numbers and Backlight 3 WHO Hypertension Risk Indicator 4 Irregular Heartbeat detector (IHB) 5 Average of the last 3 measurements 6 2 x 120 memories with date & time indication for two users | <ol style="list-style-type: none"> 7 Cuff Tightness indicator 8 Body Movement indicator 9 Low Battery indicator 10 Type-C USB port 11 22-42 cm Cuff 12 4 x "AAA" Alkaline Battery 13 TALKING in 6 languages (TR, EN, FR, ES, AR, KU) |
|--|---|



Order Code: M028739 - Product No: 6948158401161

AS9024

TENS & EMS Electrical Stimulator

- | | |
|---|--|
| <ol style="list-style-type: none"> 1 TENS – EMS – MASSAGE functions 2 2 output with 24 modes 3 20 level intensity 4 20 minutes default timer 5 Adjustable from 10-60 minutes | <ol style="list-style-type: none"> 6 Screen with backlight 7 Rechargeable lithium battery 8 Battery indicator (shows when the battery needs to be recharged) 9 Type-C charging interface |
|---|--|

Order Code: M028739 - Product No: 6970890320790



BABYGERÄUSCH

Fetal Doppler

- 1 Pocket type
- 2 Accurate Detection
- 3 Lightweight
- 4 LCD Display
- 5 Backlight
- 6 3.0 MHz Ultrasound Frequency
- 7 Battery Indicator
- 8 3 FHR modes: Manual Counting, Real-Time & Average
- 9 Adjustable Volume Level
- 10 Audible Alarm



Order Code: M028739 - Product No: 6951740520263

PROVAC

Portable Phlegm Suction Machine

Portable absorb phlegm apparatus designed based on developing orientation of similar products at home and abroad is a new generation of oil free lubrication suction device, which is suitable for use by the patient who has difficulty in phlegm removal due to illness, coma and operation, as well as for aspirating such liquid as pus and blood during the clinical practice. It is the commonly applied medical device for use in the emergency room, operation room, and for nursing in sickroom and home health care.

- 1 With a low flow, medium vacuum and fashionable stylish
- 2 Adopting oil-free piston pump
- 3 Adjustable pressure button to select any suitable pressure level
- 4 Overflow protection can safely and efficiently prevent the liquid from entering into the pump



Order Code: M028739 - Product No: 6975070911812

MOBINEB

Mesh Nebulizer

- 1 Dual Power supply with choice of Rechargeable Battery or USB Cable
- 2 Two Level Nebulization Rate Adjustment
- 3 Tiny particle size < 5 μm
- 4 Big medication cup
- 5 Low noise
- 6 Light-Weight
- 7 7 Easy to Use



Order Code: M029559 - Product No: 8698864016453



INDEX

[illegible]

DTS-3000

Digital Tourniquet System

- 1 Separate pressure and operation for 2 channels
- 2 Easy settings by touch key operation
- 3 Graphic LCD screen
- 4 Pressure and system self-testing
- 5 Self-air pressure sensor calibration
- 6 Safety alarm system (cuff, pressure, battery, time check alarms)
- 7 Built-in battery
- 8 Stand-by mode
- 9 2 channels (4 ports for 2 double or single cuffs use)
- 10 RIVA feature
- 11 LOP sensor application
- 12 Optional 14 different cuff sizes (Single, IVRA, Conic)



Order Code: M067803 - Product No: 8698864012745



LX

Electrosurgical Unit – 300 WATT

- ❶ 6SENSE Tissue Feedback Technology
- ❷ Self Control and Error Code Technology
- ❸ Wide Range of Surgical Application
- ❹ 2 Bipolar, 7 Monopolar Modes
- ❺ Endoscopic (endocut) mode
- ❻ 1 Monopolar, 1 Bipolar Output
- ❼ 10 User Adjustable Programs
- ❽ Patient Return Electrode Monitoring
- ❾ Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- ❿ Audible and Visual Alarms



Order Code: M066851 - Product No: 8698864013643

EASY +

Electrosurgical Unit – 300 WATT

- ❶ 6SENSE Tissue Feedback Technology
- ❷ Self Control and Error Code Technology
- ❸ Wide Range of Surgical Application
- ❹ 2 Bipolar, 7 Monopolar Modes
- ❺ Endoscopic (endocut) mode
- ❻ 1 Monopolar, 1 Bipolar Output
- ❼ 10 User Adjustable Programs
- ❽ Patient Return Electrode Monitoring
- ❾ Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- ❿ Audible and Visual Alarms



Order Code: M066852 - Product No: 8698864013650

MX

Electrosurgical Unit – 300 WATT

- ❶ 6SENSE Tissue Feedback Technology
- ❷ Self Control and Error Code Technology
- ❸ Wide Range of Surgical Application
- ❹ 4 Bipolar, 7 Monopolar Modes
- ❺ Endoscopic (endocut) mode
- ❻ 2 Monopolar, 1 Bipolar Output
- ❼ Simultaneous Coagulation
- ❽ 30 User Adjustable Programs
- ❾ Patient Return Electrode Monitoring
- ❿ Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- ⓫ Audible and Visual Alarms

Order Code: M066850 - Product No: 8698864013636



SIGMA +

Electrosurgical Unit – 400 WATT

- ❶ 4,7" Colour LCD Display
- ❷ 6SENSE Tissue Feedback Technology
- ❸ Self Control and Error Code Technology
- ❹ Wide Range of Surgical Application
- ❺ 4 Bipolar, 8 Monopolar Modes
- ❻ Endoscopic (endocut) mode
- ❼ 2 Monopolar, 1 Bipolar Output
- ❽ Simultaneous Coagulation
- ❾ 30 User Adjustable Programs
- ❿ Patient Return Electrode Monitoring
- ⓫ Upgradability: Beam Coagulation / Ultrasonic Aspirator/ Smoke Evacuator
- ⓬ Audible and Visual Alarms

Order Code: M066853 - Product No: 8698864013667



STAR5000

Fetal & Maternal Monitor

- ❶ 12.1" Color TFT LCD touch screen
- ❷ Advanced ARM system, convenient for maintenance and upgrade
- ❸ Standard and Big Font parameter interfaces
- ❹ High sensitivity waterproof probes
- ❺ Built-in 112mm width thermal printer
- ❻ Full data review up to 120 hours
- ❼ Automatic fetal status scoring system
- ❽ Wall-mountable
- ❾ Wired/Wireless central monitoring system options up to 128 beds
- ❿ Option of connectivity to external printer
- ⓫ 3 levels alarm system
- ⓬ Standard configurations: FHR, TOCO, FM, AFM, Thermal printer
- ⓭ Optional configurations: FHR2, Fetal Stimulator, External Printer

Order Code: M066201 - **Product No:** 8698864013629



STAR5000 E

Fetal & Maternal Monitor

- ❶ 5.7 inch TFT touchscreen
- ❷ Built-in 112mm width thermal printer
- ❸ Built-in Li-ion battery, up to 4 hours operating time
- ❹ Manual / Automatic FM measurement
- ❺ Automatic fetal status scoring system
- ❻ FHR signal strength indicator
- ❼ 60 hours real-time waveform reply & storage
- ❽ Timing monitoring function to avoiding over-time monitoring
- ❾ Optional connectivity to central monitoring system

Order Code: M066200 - **Product No:** 8698864013612



CLARMAC/CLARMIL

Standard Laryngoscope Set - Macintosh/Miller

- 1 Within plastic case
- 2 Autoclavable at 134°C
- 3 Curved and straight blade types
- 4 Resistant to disinfectant solution
- 5 The set is composed of:
 - 3 pieces of blade (Sizes: 00, 0, 1, 2, 3, 4)
 - 1 piece of handle with special lock system
 - 1 piece of spare bulb
 - Optional in sets of 3, 4, 5 blades

CLARMAC - Product No: 8698864008595

CLARMIL - Product No: 8698864008625



FLARMAC/FLARMIL

Fiber Optic Laryngoscope Set - Macintosh/Miller

- 1 Within plastic case
- 2 Autoclavable at 134°C
- 3 Curved and straight blade types
- 4 Resistant to disinfectant solution
- 5 The set is composed of:
 - 3 pieces of blade (Sizes: 00, 0, 1, 2, 3, 4)
 - 1 piece of handle with special lock system
 - 1 piece of spare bulb
 - Optional in sets of 3, 4, 5 blades

FLARMAC - Product No: 8698864008649

FLARMIL - Product No: 8698864008663



DIASET

Diagnostic Set - Standard/Halogen

- 1 Within zipper case
- 2 Autoclavable at 134°C
- 3 Adjustable light intensity
- 4 The set is composed of:
 - 1 piece of otoscope
 - 1 piece of ophthalmoscope
 - 2 pieces of larynx mirror (sizes of 3 and 4)
 - 1 piece of metal tongue depressor
 - 1 piece of adapter for disposable tongue depressor
 - 3 pieces of ear speculum (sizes of 2,5/3,5/4,5 mm)
 - 1 piece of light holder (for larynx mirrors)
 - 2 pieces of auxiliary bulb
 - 1 piece of nose speculum (reusable)
 - 1 piece of handle with screw lock system



Order Code: M040410 - **Product No:** 8698864008687

OTOSET

Otoscope Set / Ensemble D'otoscope - Standard/Halogen

- ① Within zipper pouch
- ② Adjustable light intensity
- ③ The set is composed of:
 - 1 piece of otoscope head with screw lock system
 - 3 pieces of reusable ear speculum (2,5/3,5/4,5 mm)
 - 1 piece of handle with screw lock system

Order Code: M040411 - **Product No:** 8698864008700



OFTASET

Ophtalmoscope Set / Ensemble D'ophtalmoscope - Standard/Halogen

- ① Within plastic case
- ② Autoclavable at 134°C
- ③ 24 diopter setting between -25 and +40
- ④ The set is composed of:
 - 1 piece of ophtalmoscope head with screwed lock system
 - 1 piece of handle with screw lock system

Order Code: M040412 - **Product No:** 8698864008700



OTO-OFTASET

Oto/Ophtalmoscope Set Ensemble De Oscope + Ophtalmoscope
Standard/Halogen

- ① Within plastic case
- ② Adjustable light intensity
- ③ Vacuum type bulb
- ④ 24 diopter setting at a range between -25 and +40
- ④ The set is composed of:
 - 3 pieces of reusable ear speculum (2,5/3,5/4,5 mm)
 - 1 piece of otoscope head
 - 1 piece of ophtalmoscope head
 - 1 piece of handle with screw lock system

Order Code: M040413 - **Product No:** 8698864008724





INDEX

[illegible]

pM-6000

Patient Monitor - Compact Design

- ❶ 8.4" LCD color screen
- ❷ 2 hours of battery life
- ❸ Simultaneous 8 waveforms displaying
- ❹ Visual and audible alarms
- ❺ 72 hours of trend saving
- ❻ Standard parameters: ECG, SpO2, TEMP, NIBP, HR, RESP
- ❼ Optional parameters: 2 Channels IBP, CO2

Order Code: M068101 - Product No: 8698864008380



pM-6800

Patient Monitor - Compact Design

- ❶ 8" LCD color screen
- ❷ 2 hours of battery life
- ❸ Simultaneous 7 waveforms displaying
- ❹ Visual and audible alarms
- ❺ 72 hours of trend saving
- ❻ Apnea alarm
- ❼ Standard parameters: ECG, SpO2, TEMP, NIBP, HR, RESP

Order Code: M068102 - Product No: 8698864008397



pM-8000

Patient Monitor - Compact Design

- ❶ 12.1" LCD color screen
- ❷ 2 hours of battery life
- ❸ Simultaneous 8 waveforms displaying
- ❹ Visual and audible alarms
- ❺ 72 hours of trend saving
- ❻ Apnea alarm
- ❼ Standard parameters: ECG, SpO2, TEMP, NIBP, HR, RESP
- ❼ Optional parameters: 2xIBP, ETCO2

Order Code: M068103 - Product No: 8698864008373



ECG300G

ECG Device - 3 Channels / Tabletop

- 1 Lightweight and compact design
- 2 Print modes: Automatic, manual and rhythm
- 3 Recording system: Thermal printer
- 4 AC, EMG and DFT filters
- 5 Internal battery
- 6 Internal memory: 1.000 patients
- 7 Lead placement diagram
- 8 Display: 320x240 color LCD



Order Code: M066923 - **Product No:** 6945040100089

ECG600G

ECG Device - 6 Channels/Tabletop/Touch Screen

- 1 Lightweight and compact design
- 2 Print modes: Automatic, manual and rhythm
- 3 Recording system: Thermal printer
- 4 AC, EMG and DFT filters
- 5 Internal memory: 1.000 patients
- 6 Lead placement diagram
- 7 Display: 800x480 color LCD touch screen
- 8 Interpretation function



Order Code: M066925 - **Product No:** 6945040100485

ECG1200G

ECG Device - 12 Channels/Tabletop/Touch Screen

- 1 Lightweight and compact design
- 2 Print modes: Automatic, manual and rhythm
- 3 Recording system: Thermal printer
- 4 AC, EMG and DFT filters
- 5 Internal battery
- 6 Internal memory: 1.000 patients
- 7 Lead placement diagram
- 8 Display: 800x600 color LCD touch screen



Order Code: M066926 - **Product No:** 6945040100492

EH-12

ECG Holter System - 12 Channels

- 1 Lightweight and compact design
- 2 12 channels ECG recording
- 3 1.2 GB built-in memory
- 4 Analysis channel selection
- 5 Patient event marking
- 6 Template, ST, HRV, QT, Atrial Fibrillation, HRT and T-Wave analyzes
- 7 Advanced reporting

Order Code: M067121 - **Product No:** 869886400885



ABPM-50

Ambulatory Blood Pressure Holter System - Large and Color Screen

- 1 Compact and lightweight design
- 2 Oscillometric measurement
- 3 Windows based PC software
- 4 Programmable software for day and night
- 5 Advanced reporting
- 6 Audible and visual alarms

Order Code: M067102 - **Product No:** 8698864012608



SP-10

Spirometry System - Portable

- 1 1.8" TFT color screen
- 2 Test items: FVC
- 3 Internal battery

Order Code: M066910 - **Product No:** 8698864009325



PLASMED

INDEX

SURGICAL BLADES	43
RESUSCITATORS	44
SKIN STAPLER & REMOVER	45
ELECTROSURGICAL PADS	45
INFUSION SETS	46
I.V. CANNULA	46
ECG ELECTRODE	47
LATEX EXAMINATION GLOVES	47
LATEX SURGICAL GLOVES	47

- 1 Carbon steel
- 2 Sizes: 10, 11, 12, 15, 18, 20, 21, 22, 23, 24
- 3 Individual sterile package
- 4 Packing: 100 pcs/box, 5,000 pcs/carton



NO:10



NO:11



NO:12



NO:15



NO:18



NO:20



NO:21



NO:22



NO:23



NO:24

NO:10 - Order Code: M035037 - Product No: 6945630118821

NO:11 - Order Code: M035038 - Product No: 6945630118838

NO:12 - Order Code: M035039 - Product No: 6945630118845

NO:15 - Order Code: M035040 - Product No: 6945630118852

NO:18 - Order Code: M035041 - Product No: 6945630118869

NO:20 - Order Code: M035042 - Product No: 6945630118876

NO:21 - Order Code: M035043 - Product No: 6945630118883

NO:22 - Order Code: M035044 - Product No: 6945630118890

NO:23 - Order Code: M035045 - Product No: 6945630118906

NO:24 - Order Code: M035046 - Product No: 6945630118913

AMBUSET-D

PVC Resuscitator Kit - Disposable

- 1 Pressure relieving valve
- 2 Inflatable aircushion mask with syringe
- 3 Oxygen balloon volumes: Adult (1700ml), pediatric (600ml) and infant (280ml)
- 4 Oxygen reservoir volumes: Adult (2000 ml), pediatric and infant (1600 ml)
- 5 1 piece PVC mask, 3 pieces airway and 1 piece nasal oxygen cannula in kit
- 6 Packing: 24 pcs/carton

D280 - Order Code: M034523 - Product No: 8698864013193

D600 - Order Code: M034524 - Product No: 8698864013209

D1700 - Order Code: M034525 - Product No: 8698864013216



AMBUSET-R

Silicone Resuscitator Kit - Reusable

- 1 Autoclavable (134°C)
- 2 Pressure relieving valve
- 3 Oxygen balloon volumes: Adult (1700ml), pediatric (600ml) and infant (280ml)
- 4 Oxygen reservoir volumes: Adult (2000 ml), pediatric and infant (1600 ml)
- 5 Mask sizes: 0, 1, 2, 3, 4, 5
- 6 2 pieces silicone mask, 3 pieces airway, 1 piece nasal oxygen cannula and 1 piece mouth opener in kit
- 7 Hard plastic carrying case in order to carry all the equipments
- 8 Packing: 12 pcs/carton

R280 - Order Code: M034501 - Product No: 8698864013162

R600 - Order Code: M034521 - Product No: 8698864013179

R1700 - Order Code: M034522 - Product No: 8698864013186



F-35W

Stapler Remover / Reusable / Sterile

- ❶ Can be sterilized by ethylene oxide (ETO)
- ❷ Individual sterile blister package
- ❸ Packing: 200 pcs/carton

Order Code: M035031 - **Product No:** 5060186750125



C-WR

Skin Stapler / Disposable / Sterile

- ❶ Stainless steel
- ❷ 35 pcs wide (w) hooks
- ❸ Individual sterile blister package
- ❹ Packing: 5 pcs/box, 100 pcs/carton

Order Code: M035031 - **Product No:** 5060186750125



GP03 / GP04

Electrosurgical Pad - Disposable

- ❶ Monopolar or bipolar
- ❷ Hydrogel conductive adhesive
- ❸ Compatible with all electrosurgical units
- ❹ Reusable cables (Optional)
- ❺ Packing: 50 pcs/box, 400 pcs/carton

GP03 - Order Code: M034302 - **Product No:** 8698864001404

GP04 - Order Code: M034304 - **Product No:** 8698864001428



INFSET-S

Infusion Set (Vented) Disposable / Sterile

- ❶ For infusion bottles and bags
- ❷ Vented spike
- ❸ Air-inlet filter with snap on cap
- ❹ 20 drops= 1 ml
- ❺ 15 µm filter
- ❻ Roller clamp
- ❼ injection site
- ❽ Luer slip connector
- ❾ 21 G needle
- ❿ Tube length: 150 cm
- ⓫ Individual sterile package
- ⓬ Packing: 500 pcs/carton


Order Code: M034231 - **Product No:** 8698864009912



PLASCATH

I.V. Cannula

- ❶ Japanese Stainless Steel Cannula
- ❷ Extra Fine Needle Point for Smooth Cannulation
- ❸ Thin Walled FEP/PTFE Catheter for Lower Penetration Force
- ❹ Customized Tipping Technology for Lower Tip-Wall Thickness
- ❺ Minimum Clearance between Catheter and Needle for Zero Peel-Back Issues during use
- ❻ Luer Lock Designed to Prevent Leakage
- ❼ Non-pyrogenic
- ❽ No: 14G, 16G, 17G, 18G, 20G, 22G, 24G With Wings & Injection Port Radio-opaque Catheter for Adult
- ❾ 24G, 26G With Small Wings & Without Injection Port
- ❿ Radioopaque Catheter for Neonatal



No:14G Adult - Order Code: M031825 - **Product No:** 8698864014671

No:16G Adult - Order Code: M031826 - **Product No:** 8698864014688

No:17G Adult - Order Code: M031827 - **Product No:** 8698864014695

No:18G Adult - Order Code: M031828 - **Product No:** 8698864014701

No:20G Adult - Order Code: M031829 - **Product No:** 8698864014718

No:22G Adult - Order Code: M031830 - **Product No:** 8698864014725

No:24G Adult - Order Code: M031831 - **Product No:** 8698864014732

No:24G Neonatal - Order Code: M031841 - **Product No:** 8698864014749

No:26G Neonatal - Order Code: M031842 - **Product No:** 8698864014756

ECG ELECTRODE

- 1 B5055 Adult: Round, Size: 50*55mm, Foam Backing, Solid Gel, 50 pcs / pouch.
- 2 B2525 Size: Child: 25.4*25.4mm, Foam Backing, Solid Gel, 50 pcs / pouch.

B5055 - Order Code: M035905 - Product No: 8698864016187

B2525 - Order Code: M035906 - Product No: 8698864016194



LATEX EXAMINATION GLOVES

- 1 Ambidextrous
- 2 Single use only
- 3 Powdered / Powder-free
- 4 Sizes: X-Small, Small, Medium, Large, X-Large



POWDERED / X-SMALL - Product No: 8698864013834

POWDERED / SMALL - Product No: 8698864013841

POWDERED / MEDIUM - Product No: 8698864013858

POWDERED / LARGE - Product No: 8698864013865

POWDERED / X-LARGE - Product No: 8698864013872

POWDERED-FREE / X-SMALL - Product No: 8698864013889

POWDERED-FREE / SMALL - Product No: 8698864013896

POWDERED-FREE / MEDIUM - Product No: 8698864013902

POWDERED-FREE / LARGE - Product No: 8698864013919

POWDERED-FREE / X-LARGE - Product No: 8698864013926

LATEX SURGICAL GLOVES

- 1 Natural rubber latex
- 2 Single use only
- 3 Powdered / Powder-free
- 4 Sizes : 6 - 6,5 - 7 - 7,5 - 8 - 8,5



POWDERED / NO:6 - Product No: 8698864014763

POWDERED / NO:6,5 - Product No: 8698864014770

POWDERED / NO:7 - Product No: 8698864014787

POWDERED / NO:7,5 Product No: 8698864014794

POWDERED / NO:8 - Product No: 8698864014800

POWDERED / NO:8,5 - Product No: 8698864014817

POWDERED-FREE / NO:6 - Product No: 8698864014824

POWDERED-FREE / NO:6,5 - Product No: 8698864014831

POWDERED-FREE / NO:7 - Product No: 8698864014848

POWDERED-FREE / NO:7,5 - Product No: 8698864014855

POWDERED-FREE / NO:8 - Product No: 8698864014862

POWDERED-FREE / NO:8,5 - Product No: 8698864014879

DISTRIBUTION PRODUCTS

INDEX

MP-60

Infusion Pump - Modular Design

- | | |
|------------------------------|---------------------------------|
| 1 Touchscreen | 9 Adjustable KVO |
| 2 Lightweight and small | 10 Night mode |
| 3 Modular desing | 11 Patient data |
| 4 Plug&play | 12 Drug library |
| 5 Stackable | 13 IV standına sabitleme imkânı |
| 6 Automatic door | 14 Fixing to IV pole |
| 7 Anti free flow mechanism | 15 Audible and visual alarms |
| 8 Automatic and manual bolus | 16 Internal battery |

MEDCAPTAIN



Order Code: M066830 - Product No: 8698864012639

SYS-6010

Infusion Pump – Portable

- 1 User-friendly
- 2 Touchscreen
- 3 Anti free flow mechanism
- 4 High safety and stability
- 5 High accuracy
- 6 Optional drop sensor

MEDCAPTAIN



Order Code: M066831 - Product No: 8698864012646

MP-80

Infusion Workstation

- 1 Modular desing
- 2 Stackable
- 3 Plug&Play
- 4 Synchronize alarm with pumps
- 5 Wi-Fi Connectivity to central monitoring system
- 6 Compatible with MP-60 infusion pump and MP-30 syringe pump

MEDCAPTAIN



Order Code: M066834 - Product No: 8698864012653

MP-30

Syringe Pump - Modular Design

- | | |
|------------------------------|------------------------------|
| 1 Touchscreen | 8 Night mode |
| 2 Lightweight and small | 9 Patient data |
| 3 Modular desing | 10 Drug library |
| 4 Plug&play | 11 Fixing to IV pole |
| 5 Stackable | 12 Audible and visual alarms |
| 6 Automatic and manual bolus | 13 Internal battery |
| 7 Adjustable KVO | |

Order Code: M066832 - Product No: 8698864012615

MEDCAPTAIN



SYS-3011

Syringe Pump - Portable

- 1 User-friendly
- 2 High safety and stability
- 3 High accuracy
- 4 Patented drive system
- 5 Special technology of motor drive

Order Code: M066833 - Product No: 8698864012622

MEDCAPTAIN



HI-101

Spirometry System - Table Top

- | | |
|--------------------------------------|---|
| 1 For testing asthma, COPD and other | 10 Just pres the buttons and testing is |
| 2 respiratory failure patients | completed in 3 minutes |
| 3 Compact for mobility | 11 Rechargeable lead acid battery |
| 4 Large LCD screen | 12 Thermal printer |
| 5 Test items: SVC, FVC, MVV, BD | 13 Data managment software for Windows |
| 6 Proven ATS accuracy | 14 Optional items: Calibration syringe |
| 7 Data storage: Approx. 100 tests | 15 (3L or 1 L), spiro filtre |
| 8 Computer interface: RS-232C | |
| 9 Ideal for clinics | |

Order Code: M063901 - Product No: 8698864013353

CHEST



BIOPSY PUNCH

Sterile

- 1 User-friendly
- 2 High safety and stability
- 3 High accuracy
- 4 Patented drive system
- 5 Special technology of motor drive



2 mm - Order Code: M032307 - Product No: 4560146922605
3 mm - Order Code: M032301 - Product No: 4560146922599
4 mm - Order Code: M032302 - Product No: 4560146922582
5 mm - Order Code: M032303 - Product No: 4560146922612
6 mm - Order Code: M032304 - Product No: 4560146922575

Male External (Condom) Catheter

Disposable / Sterile

- 1 Soft, pure latex
- 2 Compatibility to all standard urine bags
- 3 With soft adhesive tapes
- 4 Sizes: Medium (25 mm), Large (30 mm) and Extra Large (35 mm)
- 5 Individual sterile blister package
- 6 Packing: 50 pcs/box, 600 pcs/carton



CATHETER-L - Order Code: M030716 - Product No: 8905159006694
CATHETER-M - Order Code: M030717 - Product No: 8905159006687
CATHETER-XL - Order Code: M030710 - Product No: 8905159006700

NEW ASKIR - 230/12V BR

Battery-Operated Suction Unit

NEW ASKIR 230/12V BR is an electric suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. Suitable for emergency applications and for post-operative therapy, both in professional environment and home-care. Three different options for operation: AC/DC adapter, rechargeable battery and 12V make it a versatile suction unit approved for NON-STOP operation without overheating. Provided with visual and acoustic alarm indicating low battery level. The newly-designed ergonomic handle ensures a comfortable and secure grip.

Order Code: M066933 - Product No: 8054619100370



Various TRIMPEKS products which do not take place in the Product Catalogue

Hospital Furniture

Hospitalisation Bed
Baby Cot
Mayo Table
Nightstand
Delivery Table
Operating Table
Examination Table
Blood Sampling Chair
Dialysis Chair
Ward Screen
Couch Steps
Stretcher
Trolleys / Carts
Mobile Infusion Stand (Serum holder)
Adjustable Stool
Film Viewer
Surgical Sink
Hospital Waiting Bench

Miscellaneous

ENT Consultation Unit
Dialysis Machine
Medical Refrigerator
ICU Ventilator
Emergency & Transport Ventilator
Anesthesia Machine
Defibrillator
AED (Automated External Defibrillator)
Endoscopy Column
Laparoscopy Set
Feces Analyzer
Dental Chair

Radiology Products

Mobile X-Ray
Ultrasound
Dental Panoramic X-Ray
CT Scanner

Lamps

Examination Lamp
Surgical Lamp

Digital and Mechanical Scales

Digital Baby Scale
Mechanical Baby Scale with sliding weights
Mechanical Adult Scale with sliding weights and Height Rod

Consumables

IV Cannula
Infusion Set
Disposable Syringes
Angiographic Syringes
Plasters, Bandages, Compresses, Dressings, Gauze Rolls, Cotton, etc.

RUTF (Ready to Use Therapeutic Food)

Produits divers de TRIMPEKS qui ne figurent pas dans le Catalogue de Produits

Mobilier Hospitalier

Lit d'hospitalisation
Berceau
Table de Mayo
Table de chevet
Table d'accouchement
Table d'opération
Chaise / Lit / Table / Divan d'Examen
Fauteil de prelevement
Fauteil de dialyse
Paravent
Escabeau
Brancard
Chariots
Porte-sérum (Potence)
Tabouret réglables
Négatoscope
Lavabo chirurgical
Banc d'attente

Divers

Unité de consultation ORL
Generateur de dialyse
Réfrigérateur médicalisé
Ventilateur de soins intensifs, Ventilateur d'urgence et de transport
Appareil d'anesthésie
Défibrillateur
Appareil AED
Colonne d'endoscopie
Ensemble de coelioscopie
Analyseur fécal
Fauteil dentaire

Radiologie

Radiographie Mobile
Ecographe
Radio Panoramique
Scanner CT

Lampes

Lampe d'examen
Lampe scialytique

Balances numériques et mécaniques

Pèse-bébé numérique
Pèse bébé mécanique à poids coulissant
Pèse personne mécanique à poids coulissant

Consommables

Cathéter intraveineuse
Set de perfusion intraveineuse
Seringues jetables
Seringues angiographiques
Sparadraps, Rouloaux de Gaze, Bandages, Compresses, Pansements, Cotton, etc.,

ATPE (Aliments thérapeutiques prêts à l'emploi)

TRIMPEKS

Your
“One-Stop Shop”
for Healthcare
Requirements for

- ✓ **Detection**
- ✓ **Risk Assessment**
- ✓ **Diagnosis**
- ✓ **Therapy**



Contact Us



Sultan Selim Mah. Yunus Emre Cd. No:1/11 34415 Kağıthane İstanbul Türkiye



www.trimpeks.com



export@trimpeks.com



+90 212 319 50 00